



Anti-HCVgp1 (aa 1-77) polyclonal antibody (CABT-BL1764)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment: IPASAYQVRN ASGLYHVTND CSNSSIVYEA AGMIMHTPGC VPCVRENNAS RCWVALTPTL AARNTSIPTT TIRRHVD, corresponding to amino acids 1-77 of Hepatitis C Virus Hepatitis C Virus
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Virus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Virus surface
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [951475](#)

