



Anti-HCV gp1 polyclonal antibody (CABT-BL1765)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to N-terminal residues of Hepatitis C Virus E1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HCV
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Endoplasmic reticulum membrane; Single-pass type I membrane protein
Format	Liquid
Size	100 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [951475](#)

