



Anti-HHBVgp4 polyclonal antibody (CABT-BL1788)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Hepatitis B Virus Surface Antigen
Isotype	IgG
Source/Host	Horse
Species Reactivity	HBV
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Virion membrane
Format	Liquid
Buffer	0.01M PBS
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	2703546
Protein Refseq	NP_040999
UniProt ID	P13847