



Anti-HHEX (aa 211-260) polyclonal antibody (CABT-BL1792)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 211-260 (SLDTSCEQGQ DLPSEQNKGA SLDRSQCSPS PASQEDPDSE ISESDQEVD) of Mouse Hex.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	This gene encodes a member of the homeobox family of transcription factors, many of which are involved in developmental processes. Expression in specific hematopoietic lineages
---------------------	--

suggests that this protein may play a role in hematopoietic differentiation. [provided by RefSeq, Jul 2008]

GENE INFORMATION

Entrez Gene ID	15242
----------------	-----------------------

Protein Refseq	NP_032271
----------------	---------------------------

UniProt ID	P43120
------------	------------------------
