



Anti-HIC1 (aa 400-500) polyclonal antibody (CABT-BL1797)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 400 - 500 of Human HIC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Nucleus.
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	This gene functions as a growth regulatory and tumor repressor gene. Hypermethylation or deletion of the region of this gene have been associated with tumors and the contiguous-gene syndrome, Miller-Dieker syndrome. Alternative splicing of this gene results in multiple transcript
--------------	--

GENE INFORMATION

Entrez Gene ID [3090](#)

Protein Refseq [NP_001091672](#)

UniProt ID [Q14526](#)
