



Anti-HMGN4 (full length) polyclonal antibody (CABT-BL1819)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full length Human HMGN4 (amino acids 1-90) (BC001282).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nucleus.
Format	Liquid
Buffer	1% BSA, 99% PBS pH: 7.20
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The protein encoded by this gene, a member of the HMGN protein family, is thought to reduce the compactness of the chromatin fiber in nucleosomes, thereby enhancing transcription from chromatin templates. Transcript variants utilizing alternative polyadenylation signals exist for this gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [10473](#)

Protein Refseq [NP_006344](#)

UniProt ID [A0A024R046](#)
