



Anti-HNRNPUL2 (internal region) polyclonal antibody (CABT-BL1824)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to internal of human HNRNPUL2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	IHC, ELISA
Format	Liquid
Buffer	In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	221092
Protein Refseq	NP_001073027
UniProt ID	Q1KMD3