



# Anti-HRNR polyclonal antibody (CABT-BL1848)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human HRNR.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	HRNR (hornerin) is a protein-coding gene. Diseases associated with HRNR include irritant dermatitis, and dermatitis. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is RPTN.
--------------	--

## GENE INFORMATION

Entrez Gene ID [388697](#)

---

Protein Refseq [NP\\_001009931](#)

---

UniProt ID [Q86YZ3](#)

---