



# Anti-HSF3 (aa 171-220) polyclonal antibody (CABT-BL1859)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within internal sequence amino acids 171-220 (EDKYKSLLDR VMPILKESKK LISSGDQPSG DDGEHPKVPV QDKPMNEESL) of Mouse B230358A15Rik.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Isoform 1: Cytoplasm. Nucleus. Note: Cytoplasmic under normal conditions. Translocates to the nucleus in response to heat shock. Isoform 2: Cytoplasm. Note: Does not translocate to the nucleus in response to heat shock.
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [245525](#)

---

Protein Refseq [NP\\_766519](#)

---

UniProt ID [Q8BL67](#)

---