



Anti-HSFX1 (aa 228-273) polyclonal antibody (CABT-BL1862)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between amino acids 228-273 of human LW1. (Note: the exact amino acid sequence is proprietary).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC-P
Cellular Localization	Nucleus (Potential). Cytoplasm.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Located on chromosome X, this gene encodes hypothetical protein LOC51402.
---------------------	---

GENE INFORMATION

Entrez Gene ID

[51402](#)
