



# Anti-HSFX1 (N-terminal) polyclonal antibody (CABT-BL1863)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from the N terminal of human LW1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nucleus (Potential). Cytoplasm.
Format	Liquid
Buffer	50% Glycerol, PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction** Located on chromosome X, this gene encodes hypothetical protein LOC51402.

## GENE INFORMATION

