



Anti-HSP17.6A (full length) polyclonal antibody (CABT-BL1867)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length recombinant HSP17.6 protein produced in E. coli and purified by conventional methods (no affinity tag)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic
Format	Liquid
Size	200 µl
Buffer	Whole serum
Preservative	None
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	831076
-----------------------	------------------------

Protein Refseq [NP_196764](#)

UniProt ID [O81822](#)
