



Anti-HSP21 polyclonal antibody (CABT-BL1868)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant HSP21 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic
Format	Liquid
Size	50 µg
Buffer	Whole serum
Preservative	None
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction	chloroplast located small heat shock protein.
--------------	---

GENE INFORMATION

Entrez Gene ID [828881](#)

Protein Refseq [NP_194497](#)

UniProt ID [P31170](#)
