



Anti-HSP90C (C-terminal) polyclonal antibody (CABT-BL1870)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	C terminal 238 amino acids of Hsp90C expressed with N and C terminal hexahistidine tags in E.coli, purified with Ni-NTA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	C. reinhardtii
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID [5728553](#)

UniProt ID [Q66T67](#)