



Anti-HTR5B (C-terminal) polyclonal antibody (CABT-BL1880)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from C terminal region of Mouse 5HT5B Receptor conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Cell membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	Whole serum, 1X PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [15564](#)

Protein Refseq [NP_034613](#)

UniProt ID [G3X9C6](#)
