



Anti-IER5L (aa 331-380) polyclonal antibody (CABT-BL1886)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 331-380 (TPFAPCKRAR FEDFCPDSSP DASNISNLIS IFGSGFSGLV SRQPDSSSEQP) of Mouse IER5L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	98% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	72500
Protein Refseq	NP_084520
UniProt ID	Q99J55