



Anti-IFI204 (internal region) polyclonal antibody (CABT-BL1888)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 101-150 (EAGPATPTST TSHMLASERG ETSATQEETS TAQAGTSTAQ ARTSTAQAGT) of Mouse IFI16.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [15951](#)

Protein Refseq [NP_032355](#)

UniProt ID [E9QKL6](#)
