



Anti-IFI27L2 (aa 81-130) polyclonal antibody (CABT-BL1889)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 81-130 (SVGSVLGACL GNSPSSSLPA EPEAKEDEAR ENVPQGEPPK PPLKSEKHEE) of Human IFI27L2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Membrane; Single-pass type I membrane protein.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [83982](#)

Protein Refseq [NP_114425](#)

UniProt ID [Q9H2X8](#)
