



Anti-ify-1 (aa 47-146) polyclonal antibody (CABT-BL1896)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: SSKTPGGKSLQSLASACKVPPSTKKNTIPIAFECYEDET DQIADVATIKKTEKHPCSPIDTANRCETFDLSAADIEDDM LNLEDQDVVLSEDRPYGDVI, corresponding to amino acids 47/146 of C. elegans IFY-1
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Caenorhabditis elegans
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	173872
Protein Refseq	NP_494931

UniProt ID

[Q18235](#)
