



Anti-IGSF3 polyclonal antibody (CABT-BL1907)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human IGSF3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	IGSF3, Immunoglobulin superfamily member 3, contains 8 Ig-like C2-type (immunoglobulin-like) domains. It is expressed in a wide range of tissues with high expression in placenta, kidney and lung.
---------------------	---

GENE INFORMATION

Entrez Gene ID [3321](#)

Protein Refseq [NP_001007238](#)

UniProt ID [O75054](#)
