



Anti-IL20RA (aa 159-175) polyclonal antibody (CABT-BL1920)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: EKWKRPEDLPVSMQQL, corresponding to amino acids 159-175 of Human IL20 Receptor alpha
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB, FC, IHC-P
Cellular Localization	Membrane.
Format	Liquid
Buffer	0.2% Gelatin, PBS
Preservative	0.05% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The IL20RA/IL20RB dimer is a receptor for IL19, IL20 and IL24. The IL20RA/IL10RB dimer is a receptor for IL26.
---------------------	--

GENE INFORMATION

Entrez Gene ID [53832](#)

Protein Refseq [NP_055247](#)

UniProt ID [Q9UHF4](#)
