



Anti-IL20RA (aa 272-553) polyclonal antibody (CABT-BL1921)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment, corresponding to amino acids 272-553 of Human IL20 Receptor alpha (BC113574).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Cellular Localization	Membrane.
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.2
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction The IL20RA/IL20RB dimer is a receptor for IL19, IL20 and IL24. The IL20RA/IL10RB dimer is a

receptor for IL26.

GENE INFORMATION

Entrez Gene ID	53832
----------------	-----------------------

Protein Refseq	NP_055247
----------------	---------------------------

UniProt ID	Q9UHF4
------------	------------------------
