



Anti-IME1 (aa 97-279) polyclonal Antibody (CABT-BL1932)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein, corresponding to amino acids 97-279 of <i>S. cerevisiae</i> IME1.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	<i>S. cerevisiae</i>
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	853556
Protein Refseq	NP_012627

UniProt ID

[P21190](#)
