



Anti-IMP5 (internal region) polyclonal antibody (CABT-BL1933)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 575-624 (NPGEDTTEIV TISENEATNP EDRSDSSEGW SDAHLPNEL PFIPPGASEE) of Human IMP5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	162540
Protein Refseq	NP_787078

UniProt ID

[Q8IUH8](#)
