



Anti-INCA1 (C-terminal) polyclonal antibody (CABT-BL1934)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A 15 amino acid synthetic peptide from near the C terminus of Human INCA1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	INCA1 (inhibitor of CDK, cyclin A1 interacting protein 1) is a protein-coding gene. Diseases associated with INCA1 include testicular cancer. GO annotations related to this gene include cyclin-dependent protein serine/threonine kinase inhibitor activity and cyclin binding.
---------------------	---

GENE INFORMATION

Entrez Gene ID [388324](#)

Protein Refseq [NP_001161459](#)

UniProt ID [Q0VD86](#)
