



Anti-INTS8 (aa 200-250) polyclonal antibody (CABT-BL1939)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region between residues 200 and 250 of Human FLJ20530
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Nuclear
Format	Liquid
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	INTS8 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large;subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2;(RNU2; MIM 180690)
--------------	---

GENE INFORMATION

Entrez Gene ID	55656
----------------	-----------------------

Protein Refseq	NP_060334
----------------	---------------------------

UniProt ID	Q75QN2
------------	------------------------
