



Anti-IOC4 (aa 1-30) polyclonal antibody (CABT-BL1941)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: MSEAIFQPTDIVLAKVKGFSAWPAMIIPNE, corresponding to N terminal amino acids 1/30 of <i>Saccharomyces cerevisiae</i> IOC4
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [855061](#)

Protein Refseq [NP_013758](#)

UniProt ID [Q04213](#)
