



Anti-IQCJ (aa 101-130) polyclonal antibody (CABT-BL1943)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide between C terminal amino acids 101-130 of Human IQCJ, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	654502
Protein Refseq	NP_001036171
UniProt ID	Q1A5X6