



Anti-IRGQ (aa 20-50) polyclonal antibody (CABT-BL1946)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide corresponding to a region within N terminal amino acids 20-50 from Human IRGQ.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	126298
Protein Refseq	NP_001007562
UniProt ID	Q8WZA9