



Anti-ISLR2 polyclonal antibody (CABT-BL1949)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 181-230 (HNPFHCGCGL VWLQAWAAST RVSLPEPDSI ACASPPALQG VPVYRLPALP) of Human ISLR2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Cell Membrane
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	57611
Protein Refseq	NP_001123608

