



Anti-ISLR2 (internal region) polyclonal antibody (CABT-BL1950)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from the internal region of Mouse ISLR2 conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-Fr, IHC-P
Cellular Localization	Cell Membrane
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [320563](#)

Protein Refseq [NP_001155007](#)

UniProt ID [E9Q9N7](#)
