



Anti-ISPD (aa 407-437) polyclonal antibody (CABT-BL1951)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, between amino acids 407-437 from the C-terminal region of Human ISPD.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	This gene encodes a 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase-like protein. Mutations in this gene are the cause of Walker-Warburg syndrome. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2012]
---------------------	--

GENE INFORMATION

Entrez Gene ID [729920](#)

Protein Refseq [NP_001094887](#)

UniProt ID [A4D126](#)
