



# Anti-ITIH5 polyclonal antibody (CABT-BL1970)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human ITIH5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	May act as a tumor suppressor.
--------------	--------------------------------

## GENE INFORMATION

Entrez Gene ID	<a href="#">80760</a>
----------------	-----------------------

Protein Refseq [NP\\_001001851](#)

---

UniProt ID [G5E9D8](#)

---