



# Anti-ITPKA polyclonal antibody (CABT-BL1972)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein corresponding to amino acids of human ITPKA.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Antigen affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC
<b>Format</b>	Liquid
<b>Buffer</b>	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Regulates inositol phosphate metabolism by phosphorylation of second messenger inositol 1,4,5-trisphosphate to Ins(1,3,4,5)P <sub>4</sub> . The activity of the inositol 1,4,5-trisphosphate 3-kinase is responsible for regulating the levels of a large number of inositol polyphosphates that are important in cellular signaling. Both calcium/calmodulin and protein phosphorylation mechanisms control its activity. It is also a substrate for the cyclic AMP-dependent protein kinase, calcium/calmodulin- dependent protein kinase II, and protein kinase C in vitro.
---------------------	--

# GENE INFORMATION

Entrez Gene ID [3706](#)

---

Protein Refseq [NP\\_002211](#)

---

UniProt ID [P23677](#)

---