



# Anti-KALRN polyclonal antibody (CABT-BL1980)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human KALRN.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Huntington's disease (HD), a neurodegenerative disorder characterized by loss of striatal neurons, is caused by an expansion of a polyglutamine tract in the HD protein huntingtin. This gene encodes a protein that interacts with the huntingtin-associated protein 1, which is a huntingtin binding protein that may function in vesicle trafficking. Alternatively spliced transcript
--------------	---

## GENE INFORMATION

Entrez Gene ID [8997](#)

---

Protein Refseq [NP\\_001019831](#)

---

UniProt ID [J3QSW6](#)

---