



Anti-KANSL3 (aa 401-431) polyclonal antibody (CABT-BL1981)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide between 401-431 amino acids from the Central region of Human KIAA1310.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	55683
Protein Refseq	NP_001108488
UniProt ID	Q9P2N6