



Anti-KCNQ5L (aa 880-897) polyclonal antibody (CABT-BL2018)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: CKAGESTDALSLPHVKLK, corresponding to amino acids 880-897 of Rat KCNQ5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Multi-pass membrane protein
Format	Liquid
Size	200 µg
Buffer	PBS, 1mg/ml BSA
Preservative	0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

