



Anti-KCNS3 polyclonal antibody (CABT-BL2020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to amino acids 26-224 of Human KCNS3
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC
Cellular Localization	Cell membrane. May not reach the plasma membrane but remain in an intracellular compartment in the absence of KCNB1.
Format	Liquid
Buffer	PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Potassium channel subunit. Modulates channel activity and reduces the ion flow.
---------------------	---

GENE INFORMATION

Entrez Gene ID [3790](#)

Protein Refseq [NP_002243](#)

UniProt ID [Q9BQ31](#)
