



Anti-KCNS3 polyclonal antibody (CABT-BL2021)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human KCNS3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC, IF
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Potassium channel subunit. Modulates channel activity and reduces the ion flow.
--------------	---

GENE INFORMATION

Entrez Gene ID [3790](#)

Protein Refseq [NP_002243](#)

UniProt ID [Q9BQ31](#)
