



Anti-KEN (aa 121-220) polyclonal antibody (CABT-BL2024)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: WKT TDSAYLSSKCGGLRDRADRRKQQTSPQNLEPDQTLK YEVDSDVDES RNAADFSSAFNSNDNCESAAECERSGGHNNK EEDEDDCTHKDNKSDKDTDE, corresponding to amino acids 121/220 of Drosophila melanogaster Ken and barbie.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [37785](#)

Protein Refseq [NP_523833](#)

UniProt ID [O77459](#)
