



# Anti-KIAA0100 (internal region) polyclonal antibody (CABT-BL2025)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from internal of human K0100.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC, IF, ELISA
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C/1 year

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">9703</a>
<b>Protein Refseq</b>	<a href="#">NP_055495</a>
<b>UniProt ID</b>	<a href="#">Q14667</a>