



Anti-KIAA0182 (aa 71-120) polyclonal antibody (CABT-BL2027)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 71-120 (EEPRGSSLSS ESSPVSSPAT NHSSPASTPK RVPMGPIIVP PGGHSVPSTP) of Human GSE1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Genetic suppressor element 1 is a protein that in humans is encoded by the KIAA0182 gene.

GENE INFORMATION

Entrez Gene ID [23199](#)

Protein Refseq [NP_001127945](#)

UniProt ID [Q14687](#)
