



Anti-KIAA0494 polyclonal antibody (CABT-BL2031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Format	Lyophilized Powder
Buffer	PBS with 0.02% Sodium Azide, 1%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20°C, Avoid freeze / thaw cycles

GENE INFORMATION

Entrez Gene ID	9813
Protein Refseq	NP_055589.1
UniProt ID	O75071