



Anti-KIAA0889 (aa 97-125) polyclonal antibody (CABT-BL2034)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide between N terminal amino acids 97-125 of Human C20orf117, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	Uncharacterized protein C20orf117 is a protein that in humans is encoded by the C20orf117 gene.
--------------	---

GENE INFORMATION

Entrez Gene ID [140710](#)

Protein Refseq [NP_542194](#)

UniProt ID [O94964](#)
