



# Anti-KIAA0930 (aa 288-337) polyclonal antibody (CABT-BL2035)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within internal amino acids 288-337 (ERVTSFSTPP TPERNNRPAF FSPSLKRKVP RNRIAEMKKS HSANDSEEFF) of human C22orf9.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** Uncharacterized protein C22orf9 is a protein that in humans is encoded by the C22orf9 gene.

# GENE INFORMATION

Entrez Gene ID [23313](#)

---

Protein Refseq [NP\\_056079](#)

---

UniProt ID [Q6ICG6](#)

---