



Anti-KIAA1199 (aa 1292-1321) polyclonal antibody (CABT-BL2040)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide between amino acids 1292-1321 from the C terminal region of Human KIAA1199.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	Protein KIAA1199 is a protein that in humans is encoded by the KIAA1199 gene.
--------------	---

GENE INFORMATION

Entrez Gene ID [57214](#)

Protein Refseq [NP_061159](#)

UniProt ID [Q8WUJ3](#)
