



# Anti-KIAA1804 polyclonal antibody (CABT-BL2048)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide from the N-terminal region of human MLK4 conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, IHC-P, WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

Introduction	KIAA1804 (Mixed Lineage Kinase 4) is a Protein Coding gene. GO annotations related to this gene include protein homodimerization activity and MAP kinase kinase kinase activity. An important paralog of this gene is MAP3K9.
--------------	---

## GENE INFORMATION

Entrez Gene ID [84451](#)

---

Protein Refseq [NP\\_115811](#)

---

UniProt ID [Q5TCX8](#)

---