



# Anti-KIAA1804 polyclonal antibody (CABT-BL2049)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from internal of human MAP3KL4.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C/1 year

## BACKGROUND

<b>Introduction</b>	KIAA1804 (Mixed Lineage Kinase 4) is a Protein Coding gene. GO annotations related to this gene include protein homodimerization activity and MAP kinase kinase kinase activity. An important paralog of this gene is MAP3K9.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [84451](#)

---

Protein Refseq [NP\\_115811](#)

---

UniProt ID [Q5TCX8](#)

---