



Anti-KIF2B polyclonal antibody (CABT-BL2055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 216-265 (VCVRKRPLNQ RETTLKDLDI ITVPSDNVVM VHESKQKVDL TRYLQNQTFC) of Human KIF2B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm; cytoskeleton; centrosome. Cytoplasm; cytoskeleton; spindle. Kinetochore. Note: Association with kinetochore is transient.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	84643
Protein Refseq	NP_115948

UniProt ID

[Q8N4N8](#)
