



Anti-KIF2C polyclonal antibody (CABT-BL2061)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide of human MCAK
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Preservative	None
Storage	-20°C

BACKGROUND

Introduction	<p>The protein encoded by this gene is a member of kinesin-like protein family. Proteins of this family are microtubule-dependent molecular motors that transport organelles within cells and move chromosomes during cell division. This protein is important for anaphase chromosome segregation and may be required to coordinate the onset of sister centromere separation.</p> <p>[provided by RefSeq, Jul 2008]</p>
--------------	---

GENE INFORMATION

Entrez Gene ID	11004
----------------	-----------------------

Protein Refseq [NP_006836](#)

UniProt ID [B7Z6Q6](#)
